



**Front-Line Device
for Your Infusion
Therapy Needs**



Features and Benefits

Easy-to-Use:

- ▶ Smooth and simple, single-step safety activation
- ▶ Clear plastic base provides excellent site visibility
- ▶ Dual port adaptability
- ▶ Easy, flexible placement and access

Designed for Patient Comfort:

- ▶ Smallest footprint of all BARD® Safety Infusion Sets
- ▶ Non-absorbent patient comfort pad
- ▶ Better dressing profile
- ▶ Lower insertion and penetration force compared to other non-siliconized port access needles

Safety:

- ▶ Robust safety mechanisms provide a front-line defense against needle sticks

Power Injection:

- ▶ Receive CECT scans
- ▶ Receive IV therapies



Product Code	Code Description WITHOUT Y-Site	Maximum Flow Rate
84446001	19 Ga x 0.75" (19mm)	5 ml/s
84446002	19 Ga x 1" (25mm)	5 ml/s
84446003	19 Ga x 1.5" (38mm)	5 ml/s
84446005	20 Ga x 0.75" (19mm)	5 ml/s
84446006	20 Ga x 1" (25mm)	5 ml/s
84446007	20 Ga x 1.5" (38mm)	5 ml/s
84446009	22 Ga x 0.75" (19mm)	2 ml/s
84446010	22 Ga x 1" (25mm)	2 ml/s
84446011	22 Ga x 1.5" (38mm)	2 ml/s


Product Code	Code Description WITH Y-Site	Maximum Flow Rate
84446013	19 Ga x 0.75" (19mm)	5 ml/s
84446014	19 Ga x 1" (25mm)	5 ml/s
84446015	19 Ga x 1.5" (38mm)	5 ml/s
84446017	20 Ga x 0.75" (19mm)	5 ml/s
84446018	20 Ga x 1" (25mm)	5 ml/s
84446019	20 Ga x 1.5" (38mm)	5 ml/s
84446021	22 Ga x 0.75" (19mm)	2 ml/s
84446022	22 Ga x 1" (25mm)	2 ml/s
84446023	22 Ga x 1.5" (38mm)	2 ml/s

Priming Volume for POWERLOC® Max SIS

19, 20, and 22 Ga.	with Y-Site - 0.4 ml
19, 20, and 22 Ga.	without Y-Site - 0.3 ml

Minimum Gravity Flow Rate for POWERLOC® Max SIS

19 Ga.	1680 ml/hr
20 Ga.	960 ml/hr
22 Ga.	312 ml/hr



Learn More

For more information on Safety Infusion Sets, please contact your Bard Access Systems representative. Or visit us on the web at www.bardaccess.com

Bard France S.A.S. (France, North Africa) FRANCE
Tel: + 33 1 39 30 58 58 Fax: + 33 1 39 30 58 59
service.commande@crbard.com www.bardfrance.fr

Bard de España S.A. SPAIN
Tel: + 34 93 253 7800 Fax: + 34 93 253 7834
spain.customerservice@crbard.com

C. R. Bard Portugal L.D.A. PORTUGAL
Tel: + 351 21 319 0330 Fax: + 351 21 319 0339
Bard.Portugal@crbard.com

Bard S.p.A. (Italy, Eastern Mediterranean, Middle East) ITALY
Tel: + 39 06 524 93 249 Fax: + 39 06 529 58 52
customerservice.italia@crbard.com www.barditalia.it

Bard Hellas SA GREECE
Tel: + 30 210 96 90 780 1 Fax: + 30 210 96 28 810
Bard.Hellas@crbard.com

Bard Medica S.A. SWITZERLAND
Tel: + 41 44 722 5360 Fax: + 41 44 722 5370
info-bce.access@crbard.com

Bard Czech Republic S.R.O. CZECH REPUBLIC
Tel: + 420 242 408 630 Fax: + 420 242 410 185
info-bce.access@crbard.com

C. R. BARD GmbH GERMANY
Tel: + 49 721 9445 124 Fax: + 49 721 9445 100
info-bce.access@crbard.com

Bard Poland Sp. z o.o. POLAND
Tel: + 48 22 321 0930 Toll Free: 0800 494900
Fax: + 48 22 321 0938 Toll Free: 0800 494901
info-bce.access@crbard.com

Bard Medica S.A. AUSTRIA
Tel: + 43 1 49 49 130 Fax: + 43 1 49 49 130 30
info-bce.access@crbard.com

Bard Limited UNITED KINGDOM
Tel: + 44 1293 529555 Fax: + 44 1293 552428
customer.services@crbard.com www.barduk.com

IRELAND
Tel: + 353 1 2130706 Fax: + 353 1 2130604
customer.services@crbard.com

Bard Norden AB SWEDEN
Tel: + 46 42 38 60 00 Fax: + 46 42 38 60 10
support.nordic@crbard.com

Bard Norway AS NORWAY
Tel: + 47 815 00 266 Fax: + 47 815 00 299
support.nordic@crbard.com

Bard Denmark AB DENMARK
Tel: + 45 49 26 20 08 Fax: + 45 49 26 20 09
support.nordic@crbard.com

Bard Finland Oy FINLAND
Tel: + 358 9 347 30 53 Fax: + 358 9 347 10 07
support.nordic@crbard.com

Bard Medical SA (Pty) Ltd. SOUTH AFRICA
Tel: + 27 11 524 9900 Fax: + 27 86 537 7250
www.bard.co.za

Bard Benelux N.V. BELGIUM
Tel: + 32 14 28 69 52 Fax: + 32 14 28 69 66
custservicebelgium@crbard.com

THE NETHERLANDS
Tel: + 31 30 60 00 572 Fax: + 31 30 60 22 345
custservicenederland@crbard.com